


Resuscitation Unit 104 E type



Specifications

| Atom's code. 60982 Resuscitation Unit 104 E type (Stand-alone type) | |
|--|--|
| General | Pressure indication: -2~8kPa (-20~80cmH ₂ O) |
| Suction unit | Suction pressure control range: 0~-26.7kPa (0~-200mmHg) |
| Oxygen blender/ flow meter unit | Oxygen concentration range: 21~100% Flow rate range: 0~15L/min |
| Resusci Flow unit | Maximum pressure limit (P _{MAX}) setting range: 2.0~5.9kPa (20~60cmH ₂ O) Default setting: 3.9kPa (40cmH ₂ O) |
| Accessories | Patient Circuit1 Test Bag, 50mL1 Piping Connecting Hose for Air ...1 Piping Connecting Hose for Oxygen ...1 Outlet Connecting Hose for Air/Oxygen ...1 Oxygen Supply Hose1 Face Mask for Neonatal1 Face Mask for Infant1 Oxygen Humidifier Jar Holder1 Humidifier Bottle OX-2301 Outlet Connecting Hose for Suction1 Suction Jar Holder1 Suction Canister, 500mL1 Dispo Liner, 500mL1 Suction Hose1 Silencer (for E-type)1 Operation Manual1 |

<Note> T-Piece Resuscitation Unit with Blender and Suction Unit of ejector type to make suction operation possible where no suction piping is available.

Options & Supplies

Atom's code. 61616 **MF-Rail** For F-Rail (Left & Right)

Atom's code. 61619 **Suction Canister, 500mL** For S401D

Atom's code. 60355 **Dispo Liner, 500mL**
For S401D, Suction Canister (30pcs/box)

① Atom's code. 60386 **Patient Circuit** For Resusci Flow (5pcs/box)

② Atom's code. 60298 **Test Bag, 50mL** (5pcs/pkg)

③ Atom's code. 60146 **Face Mask** For Infant (20pcs/box)

④ Atom's code. 60145 **Face Mask** For Neonatal (20pcs/box)



Atom Medical products are produced under strict production control in Japan.
We provide value and trust of the highest quality for the tiny baby's lives in the world.

ATOM MEDICAL CORPORATION

<https://www.atomed-global.com>